

# Acry Rock

ACRY ROCK teeth comply with  
ISO 3336: 1996 type 1 - type 2

## MOULD CHART



**Ruthinium®  
Group**

**Dental Manufacturing Spa**  
Via Ca' Mignola Nuova, 1699 -  
45021 Badia Polesine (RO)  
Tel. 042.551.628 042.559.05.87 - Fax  
Info@ruthinium.it  
Www.ruthinium.it



# Acry Rock

ACRY ROCK teeth comply with ISO 3336: 1996 type 1 - type 2



## SHADE GUIDE

A/1	WHITE	A 1
B/1		A 2
C/1	PINK	A 3
D/1		A 3,5
		A 4
E/1		B 1
F/1	YELLOW	B 2
F/2		B 3
G/1		B 4
H/1	BROWN	C 1
L/1		C 2
		C 3
		C 4
N/1		D 2
N/2	GREY	D 3
N/3		D 4

**S 13**

- (1) 42,3
- (2) 7,4
- (3) 8,9
- (4) 8,4

**S 14**

- (1) 42,6
- (2) 7,7
- (3) 9,2
- (4) 7,9

**S 15**

- (1) 41,0
- (2) 8,6
- (3) 10,3
- (4) 7,7

**S 17**

- (1) 44,8
- (2) 9,7
- (3) 12,0
- (4) 8,0

**S 18**

- (1) 45,2
- (2) 9,0
- (3) 11,0
- (4) 8,7

**S 25**

- (1) 43,8
- (2) 11,2
- (3) 11,2
- (4) 8,3

**S 26**

- (1) 45,7
- (2) 11,9
- (3) 11,7
- (4) 8,4

**S 32**

- (1) 42,8
- (2) 9,0
- (3) 10,7
- (4) 7,8

**S 33**

- (1) 46,5
- (2) 10,0
- (3) 11,5
- (4) 8,6

**S 37**

- (1) 42,7
- (2) 9,3
- (3) 11,3
- (4) 7,6

**S 38**

- (1) 41,4
- (2) 9,2
- (3) 11,1
- (4) 7,8

**S 39**

- (1) 43,3
- (2) 9,1
- (3) 11,2
- (4) 8,0

**S 40**

- (1) 41,4
- (2) 9,6
- (3) 11,6
- (4) 7,5

**S 41**

- (1) 44,0
- (2) 9,3
- (3) 11,4
- (4) 8,1

**S 42**

- (1) 43,4
- (2) 9,0
- (3) 11,2
- (4) 8,0

**S 47**

- (1) 41,8
- (2) 10,0
- (3) 11,6
- (4) 7,7

**S 48**

- (1) 44,7
- (2) 11,0
- (3) 12,6
- (4) 8,2

ACRY ROCK teeth comply with ISO 3336: 1996 type 1 - type 2



**TABELLA DI ARTICOLAZIONE**

ARTICULATION COMBINATION CHART

ARTIKULATIONS-TABELLE

TABLEAU D'ARTICULATION

TABLA DE ARTICULACION

OVERSIGHT OVER ARTIKULATION

TABELA DE ARTICULAÇÃO

TABELL ÖVER LEDRÖRELSE

ΠΙΝΑΚΑΣ ΑΡΘΡΩΣΕΩΝ

ARTICULATIETABEL

PURENTATAULUKKO

S 13	I 13	D 35
S 14	I 13	D 35
S 15	I 38	D 37 · D 43
S 17	I 41	D 39 · D 44
S 18	I 18	D 39 · D 44
S 25	I 18	D 39
S 26	I 46	D 39 · D 45
S 32	I 41 · I 42	D 39 · D 44
S 33	I 46	D 45
S 37	I 37	D 36 · D 42
S 38	I 38	D 37 · D 43
S 39	I 51	D 35
S 40	I 40	D 36 · D 42
S 41	I 41	D 39 · D 44
S 42	I 42	D 39 · D 44
S 47	I 18 · I 53	D 39
S 48	I 46	D 45
S 51	I 51	D 35
S 52	I 52	D 36 · D 42
S 53	I 53	D 39 · D 44
S 54	I 53	D 39 · D 44
S 55	I 66	D 45 · D 46
S 56	I 46	D 45
S 60	I 60	D 33
S 61	I 61	D 36 · D 42
S 62	I 61	D 36 · D 42
S 63	I 63	D 39 · D 44
S 64	I 64	D 41
S 65	I 65	D 41
S 66	I 66	D 44
S 67	I 67	D 45
S 68	I 68	D 46

# Acry Rock

ACRY ROCK teeth comply with ISO 3336: 1996 type 1 - type 2




<b>I 13</b>		
(1) 33,3 (2) 7,4 (3) 9,2 (4) 5,0		
<b>I 18</b>		
(1) 36,3 (2) 8,6 (3) 10,4 (4) 5,3		
<b>I 37</b>		
(1) 32,2 (2) 7,1 (3) 8,6 (4) 4,7		
<b>I 38</b>		
(1) 31,7 (2) 8,0 (3) 10,0 (4) 4,7		
<b>I 40</b>		
(1) 31,5 (2) 9,2 (3) 11,0 (4) 4,7		
<b>I 41</b>		
(1) 34,2 (2) 8,5 (3) 10,2 (4) 5,0		
<b>I 42</b>		
(1) 35,0 (2) 8,2 (3) 9,9 (4) 5,4		
<b>I 46</b>		
(1) 35,1 (2) 10,2 (3) 10,4 (4) 5,0		
<b>I 51</b>		
(1) 31,2 (2) 7,3 (3) 9,3 (4) 5,5		
<b>I 52</b>		
(1) 32,7 (2) 8,4 (3) 10,2 (4) 4,9		

<b>I 53</b>		
(1) 34,1 (2) 8,7 (3) 11,3 (4) 5,2		
<b>I 60</b>		
(1) 30,7 (2) 8,0 (3) 9,8 (4) 4,5		
<b>I 61</b>		
(1) 30,7 (2) 8,3 (3) 10,0 (4) 4,6		
<b>I 63</b>		
(1) 34,8 (2) 9,4 (3) 11,0 (4) 4,9		
<b>I 64</b>		
(1) 31,5 (2) 8,7 (3) 10,5 (4) 4,6		
<b>I 65</b>		
(1) 34,7 (2) 10,1 (3) 12,2 (4) 5,0		
<b>I 66</b>		
(1) 35,5 (2) 10,4 (3) 12,8 (4) 5,1		
<b>I 67</b>		
(1) 40,7 (2) 10,3 (3) 12,7 (4) 5,7		
<b>I 68</b>		
(1) 39,2 (2) 11,2 (3) 13,4 (4) 5,7		

ACRY ROCK teeth comply with ISO 3336: 1996 type 1 - type 2



**D33** (1) 26,9 - (2) 30,5 - (3) 7,2 - (4) 8,1



**D35** (1) 30,7 - (2) 33,0 - (3) 8,5 - (4) 5,5



**D42** (1) 28,7 - (2) 30,6 - (3) 8,2 - (4) 7,8



**D36** (1) 28,8 - (2) 31,1 - (3) 8,2 - (4) 8,6



**D43** (1) 29,6 - (2) 31,4 - (3) 8,3 - (4) 8,1



**D37** (1) 30,7 - (2) 33,2 - (3) 8,4 - (4) 8,8



**D44** (1) 31,1 - (2) 33,0 - (3) 8,5 - (4) 8,6



**D39** (1) 30,1 - (2) 35,0 - (3) 8,5 - (4) 8,5



**D45** (1) 31,0 - (2) 33,3 - (3) 8,8 - (4) 8,6



**D41** (1) 30,0 - (2) 33,5 - (3) 8,6 - (4) 8,9



**D46** (1) 33,6 - (2) 35,6 - (3) 9,6 - (4) 9,2

